


# USE OF THE RESPICLICK®



Drug : \_\_\_\_\_  
 Indication (maintenance or rescue or maintenance+rescue): \_\_\_\_\_  
 Number of inhalations: \_\_\_\_\_ Number of times per day: \_\_\_\_\_  
 Other information: \_\_\_\_\_



		Mastered	Not mastered
1	Remove the device from the foil package at the time of use: a. The device is ready to use		
2	With the colored cap closed, place the device upright		
3	Open the cap completely until it clicks a. The dose is ready to be inhaled		
4	Keep the device in an upright position		
5	Breathe out normally next to the device, emptying the lungs without forcing		
6	Place the mouthpiece in the mouth (between the lips and teeth)		
7	Breathe in through the mouth quickly and deeply: a. Be careful not to block the vent		
8	Move the device away from the mouth		
9	Respect a pause: hold your breath for 5 to 10 seconds		
10	Breathe out through the nose		
11	Close the cap in its original position		

## REMARKS

To limit the adverse effects of the medication, rinse your mouth with water (do not swallow the water).

**MAINTENANCE:**

**HOW TO KNOW IF IT IS EMPTY:**



**RQESR**

RÉSEAU QUÉBÉCOIS D'ÉDUCATION  
EN SANTÉ RESPIRATOIRE

La référence des professionnels